## ALKODEST Bieranalyse

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## ALKODEST Beer analysis

### **Pre-treatment of the beer sample**

The carbonic acid needs to be removed from this beer sample. To do this pour appr. 130ml of beer into a container which can hold 2 to 3 times this volume and can be closed (p.e. plastic bottle with screw cap). Close the container and shake vigorously in order to form foam. Open carefully and release the carbonic acid. Close again and repeat this several times until no pressure is formed after shaking. After this filter the beer through a folded filtering paper. Collect the filtered beer into another sealable container. Best is to place this container into a water bath at 20°C (eventually place the lead ring onto it). The warmer the sample the quicker the carbonic acid will be released.

### **Determining the density**

Pour minimum 100ml decarbonised beer (between 17 and 24°C) into a dry glass cylinder (or a glass cylinder rinsed with beer). Immerse a dry and clean densimeter into the glass cylinder and read out the density. Read out immediately the temperature of the sample with the thermometer inside the densimeter.

### **Determining the alcohol content**

- Fill the 100ml measuring flask with decarbonated beer. In order to have very precise readings take care that the flask is dry. If not rinse it at first with some ml of beer.
- Take the cooler unit off the distillation flask.
- Pour the 100ml into the distillation flask. Rinse the measuring flask twice with a small quantity of distilled water and add the rinse water to the liquid that is contained in the distillation flask.
- Add, if necessary, some drops of an antifoaming agent.
- Put 2 or 3 small boiling stones (reusable) in the distillation flask.
- Apply a very thin layer of joint grease onto the ground joint of the cooler unit and press it gently into the ground joint of the distillation flask.
- Rinse out the 100 ml measuring flask with distilled water before you use the plastic pipette to fill in approx. 5 ml of distilled water (= approx. two times the full capacity of the plastic pipette) and finally position the measuring flask under the end of the distillate delivery hose of the cooler unit. Increase the weight of the measuring flask with the lead ring and place the flask in a cold water bath.
- Turn on the hot plate or any other source of heat to be used and select a temperature setting that ensures that the liquid that is contained in the distillation flask will only be heated up to the temperature at which it just starts to turn to steam, so that the steam that is generated is transferred into the cooler unit. Make sure that the liquid to be analysed can in no case boil over into the measuring flask. Proceed with the distillation process until approx. 4/5 of the liquid volume originally contained in the distillation flask are distilled over into the measuring flask. The measuring flask is now filled to a level that leaves at least one fingerbreadth of empty space between the surface of the distillate and the base of the neck of the measuring flask.
- Remove the measuring flask immediately from the cooler unit and put in the stopper.
- Mix the contents of the measuring flask carefully (to be seen from the clearing of the distillate) by moving the measuring flask with the stopper gently in circles and then use the plastic pipette to re-fill the measuring flask precisely up to the level of its 100 ml adjustment mark with distilled water that has previously been brought to the required temperature
- Pour the contents of the measuring flask into the vertically adjusted dry glass cylinder. Then immerse a dry and clean alcoholmeter into that glass cylinder.
- Take the reading for the apparent alcohol content of the distillate from the built-in scale of the alcoholmeter. Take the temperature reading from the built-in thermometric scale of the alcoholmeter immediately after the reading for the apparent alcohol content.

---

### **Brouwland**

Korspelsesteenweg 86 • B-3581 Beverlo - Belgium  
Tel. +32-(0)11-40.14.08 • Fax. +32-(0)11-34.73.59  
sales@brouwland.com • www.brouwland.com