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Vinmetrica

Free SO₂ measurement in white, rose and red wines

The great thing about this kit is that you can perform a correct SO₂ analysis on white, rose and red wines.

Required materials:

- burette 25 ml
- 1 volumetric pipette with one mark 25 ml
- 1 measuring beaker, 100 ml or larger
- SO₂ reagents set
 - o SO₂ titrant
 - o acid solution
 - o reactant solution
- small dropper pipettes

Method:

- Fill the burette with SO₂ titrant and set to 0.
- Use a volumetric pipette to transfer 25 ml of wine to the measuring beaker (without air entrapment).
- Use the dropper pipette to add 2 ml of acid solution and 2 ml of reactant solution; mark the pipettes for each reagent to avoid contamination.
- Activate the measurement kit.
- Place the probe in the measuring beaker.
- Titrate (drop-by-drop) the SO₂ titrant from the burette until the measurement kit gives a signal. Gently swirl the measuring beaker during titration.
- Enter the number of ml SO₂ titrant in the formula below.

The free SO₂ is calculated using the following formula

$$\text{ppm SO}_2 = \frac{64 * V * N * 1000}{2 * S}$$

V = ml of titrant used to reach the endpoint; N = normality of the titrant. S = ml of wine sample.³ If you use 25 ml of wine sample as indicated and the titrant is 0.0156 N as in the set, then the calculation is simple.

$$\text{ppm SO}_2 = 20 * V \quad (\text{i.e. } 20 \times V)$$

Comment: this does not take into consideration reductones (such as ascorbic acid) that can be present in the wine. If you know that this is the case, then you can subtract 5 mg from the measured quantity. Please note: this is purely indicative.

³ $\frac{64 [\text{mg SO}_2/\text{mmol SO}_2] * V [\text{ml}] * N [\text{meq/ml}] * 1000 [\text{ml/l}]}{2 [\text{meq/mmol SO}_2] * S [\text{ml}]}$

013.002.SET Vinmetrica SO₂ analyser + reagents, complete set
013.002.LUX Vinmetrica SO₂ analyser, complete Deluxe set

Total SO₂ in white and red wines

Required materials:

- burette 25 ml
- 1 volumetric pipette with one mark 25 ml
- 1 measuring beaker 50 ml
- 1 measuring beaker 100 ml or larger
- 10 ml NaOH solution (do not pipette by mouth = DANGEROUS)
- SO₂ reagents set
- small dropper pipettes

Method:

- Fill the burette with SO₂ titrant and set to 0.
- Use a volumetric pipette to transfer 25 ml of wine to the measuring beaker (without air entrapment).
- Add 10 ml NaOH and allow to rest for 10 minutes. Use a calibrated measuring glass or syringe. (Do not pipette by mouth = DANGEROUS.)
- Use the dropper pipette to add 8 ml of acid solution and 2 ml of reactant solution; mark the pipettes for each reagent to avoid contamination. Gently swirl the measuring beaker during titration.
- Activate the measurement kit.
- Place the probe in the measuring beaker.
- Titrate (drop-by-drop) the SO₂ titrant from the burette until the measurement kit gives a signal.
- Enter the number of ml SO₂ titrant in the formula below.
- Use the same formula as provided above to calculate the total SO₂.

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