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**Procedure****Using the Ripper method (colour change):**

You will need: Vinotest reagents: iodic-iodat and starch-acid

Place 25 ml white wine in a graduated beaker.

Add 2 ml Glyoxal and let it sit for 5 minutes.

Add 5 to 10 ml of starch-acid indicator.

Titrate with iodic-iodat solution until the colour changes to blue.

**Formula**

Total titrated iodic-iodat in ml x 5 = reductones in mg/l

**Conclusion**

The actual free SO<sub>2</sub> quantity is equal to the measured free SO<sub>2</sub> content in mg/l minus the reductone content in mg/l.

**Procedure****Using the Vinmetrica analysis method:**

You will need: complete Vinmetrica kit

Place 25 ml white, rosé or red wine in a graduated beaker.

Add 2 ml Glyoxal and let it sit for 5 minutes.

Add 2 ml of acid solution and 2 ml of reactant solution.

Switch on the meter and put the electrode in the beaker.

Titrate with SO<sub>2</sub> titration solution until the meter emits a persistent audible signal.

**Formula**

Total titrated SO<sub>2</sub> solution in ml x 20 = reductones in mg/l

**Conclusion**

The actual free SO<sub>2</sub> quantity is equal to the measured free SO<sub>2</sub> content in mg/l minus the reductone content in mg/l.

**Attention**

If this product is not stored in optimal conditions, it has the tendency to form flakes. These flakes can be filtered out or dissolved when the product is heated at 50 to 60°C. The formation of flakes does not affect the quality.

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