Waterdichte pH-meter pHSCAN30

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Water-resistant pH meter pHSCAN30

Congratulations on your purchase. Please read the following instructions carefully before using this measuring instrument. Should you have any questions or remarks, do not hesitate to contact us.



INTRODUCTION

When opening the packaging, a small crystallisation may be visible on the electrode of the pH meter. This is normal. Please clean with demineralised water before first use.

The paraffin around the device serves to prevent leakage of the storage solution. Please remove the paraffin with warm water prior to first use.

Inserting batteries

1. Prior to use, unscrew the lid of the battery compartment.
2. Insert new batteries with correct polarity (- and +).

Keys

- **MEAS/HOLD functions:**

1. Switch the meter ON/OFF.
2. When pressed during measurement: records the pH value.
3. When pressed during calibration: exits calibration mode.

- **ENTER/°C functions:**


1. When pressed during calibration: confirms the value.
2. When pressed during measurement: reads the temperature.

- **CAL function:**

1. Press to activate the calibration mode.
2. Press to select a calibration point during calibration mode.

Before you start:

- Immerse the probe in a storage solution or tapped water for at least 30 minutes (no distilled or demineralised water).
- Then rinse with distilled or demineralised water.

USE

Automatic shutdown

If you do not press any key for 8 minutes, the meter will automatically turn off.

Calibration

Calibrate the meter before first use or after replacement of the electrode. We also advice to calibrate the meter twice a week if you frequently use the meter. If you only use the meter occasionally you need to calibrate the meter before each use.

1. Press MEAS to switch the meter on.
2. Thoroughly rinse the meter with distilled water (by stirring a couple of times in a glass). Briefly dry with filter paper.

3. Press and hold CAL until the display shows 'CAL 4-7'. Press 'CAL' again to calibrate at pH 7.00 and pH 10.00, display will show 'CAL 10-7'.

4. Press ENTER, 'CAL 7' appears.

5. Immerse the electrode approximately 2-3 cm in a pH 7.00 calibration solution and gently stir.

6. Press ENTER again. The meter starts to calibrate. The read-out will flash 3 times when it has stabilised. This is the pH value at the present temperature.

7. The calibration at pH 7.00 is completed. Subsequently 'CAL 4' appears. Repeat steps 5 and 6 as explained above, but this time with pH 4.00 buffer solution. This pH meter has a built-in function for automatic temperature compensation starting from 0 to 60°C. The pH meter will automatically compensate the measured value.


Brouwland


Attention:

- Make sure that both calibration solution and electrode are in good condition.
- Rinse and dry the probe before immersing it into another solution.
- Stirring gently with the electrode will improve the accuracy of the reading.
- A temporary error message 'Err' can appear when the measured value is out of the calibration range (e.g. if you immerse the meter in a pH 7.00 solution while you want to carry out a pH 4.00 calibration, the calibration will fail).
- Press MEAS to exit the calibration mode.

pH measurement

1. Press MEAS to turn the meter on.
2. Thoroughly rinse the meter with distilled water (by stirring a couple of times in a glass). Briefly dry with filter paper.
3. Immerse the probe into the solution and stir gently. Wait 3-5 seconds to stabilise the read-out. The measurement is finished.

HOLD-function

This function allows you to 'freeze' a certain value for later observation.

1. Press HOLD to fix the value. The display shows 'HOLD' along with the value.
2. Press HOLD again to exit this function.

Temperature measurement

Press ENTER/°C during measuring to read the temperature. Press again to return to the measuring mode.

Attention: the temperature of this meter lies between 0 and 60°C.

Note

This pH meter can only measure the pH values of liquids or semi-liquids. Never immerse the probe directly into highly acidic or highly alkaline liquids. Additionally, never measure liquids above 80°C.


Replacement of the electrode (item n° 013.114.4)

You can replace the module of the electrode. Should the meter fail to calibrate or gives fluctuating values, you need to replace the electrode.

1. Unscrew the electrode from the meter by twisting the ribbed collar anticlockwise.
2. Pull the old module from the meter.
3. Align the new module, using the gaps as a guide, so it will fit into the meter.
4. Gently push the module into place. Place and screw the ribbed collar back into place clockwise. It is necessary to calibrate the meter after replacing the electrode.

Brouwland

Korspelsesteenweg 86 • B-3581 Beverlo - Belgium
 Tél. +32 11 40 14 08 • Fax +32 11 34 73 59
 sales@brouwland.com • www.brouwland.com

SPECIFICATIONS

pH and temperature

°C	pH 4.01	pH 7.00	pH 10.01
0°C	4.01	7.12	10.33
5°C	4.01	7.09	10.25
10°C	4.00	7.06	10.18
15°C	4.00	7.04	10.11
20°C	4.00	7.02	10.05
25°C	4.01	7.00	10.01
30°C	4.01	6.99	9.95
35°C	4.02	6.98	9.92
40°C	4.03	6.97	9.88
45°C	4.04	6.97	9.85
50°C	4.06	6.97	9.82
55°C	4.07	6.98	9.80
60°C	4.09	6.98	9.77
70°C	4.12	6.99	9.73
80°C	4.16	7.00	9.69
90°C	4.20	7.02	9.66
95°C	4.22		8.84

Specifications

- pH range: pH 0.00 to 14.00
- pH accuracy: \pm pH 0.05
- Temperature range: 0 - 60°C
- Temperature accuracy: \pm 1°C
- Temperature compensation: 0 - 60°C
- Power: 3 x 1.5 V G13A batteries (not included)
- Dimensions: 180 mm x diameter 35 mm
- Weight: 200 g

APPENDIX

A pH meter is a precision measuring device that requires proper use and good maintenance.

Basic principles of a pH meter:

The basic components of a pH meter are a reference electrode and an indicator electrode. The reference electrode has a set value and thus provides a fixed reference against the indicator electrode. The indicator electrode can be easily distinguished from the reference electrode. The indicator electrode is very often bulb-shaped. In certain other models it can also be a window. This little glass bulb is very vulnerable and needs to be protected properly against all kinds of influences. May only come into contact with liquids. Therefore avoid to use fingers, tissues, paper, brushes, etc. The indicator electrode is positioned in an electrolytic environment.

The disadvantage of the pH meter is that the electrodes are unfortunately not permanent. Many electrodes are spent after only one year. The lifespan is determined not only by the frequency of use but also by the way in which they are used and maintained.

The better you protect your meter, the longer its longevity.

A few pointers:

- Rinse the electrode in demineralised water before and after each use. Make certain that no residues adhere to the electrode.
- Always keep the electrode in some storage solution.
- Never store the electrode in demineralised water.
- If crystals form (after a long storage without use), rinse with demineralised water and let it soak for one day in storage solution.
- If the measured value appears slowly, this indicates reduced electrode activity.
- If you have a refillable electrode that you wish to refill, first drain the electrode. The electrolyte is no longer fit for purpose and will always have to be replaced entirely. After refilling, allow the meter to stabilise for a few hours and then recalibrate.
- A temperature above 50°C is not good for the electrode and reduces its longevity.
- The reason for the repeated recalibrations of a pH meter is the result of the constantly reduced activity of the electrode.
- Always start calibration with a pH 7 buffer fluid. Follow with pH 4. When changing the buffer fluid, always rinse thoroughly with demineralised water.
 - Calibrate weekly in case of daily use.
 - Calibrate monthly in case of weekly use.
 - Calibrate each time in case of sporadic use.
 - Store the buffer fluids preferably in a refrigerator and allow them to reach room temperature before use.
- Do not reuse used buffer fluid; this will be contaminated by exposure to air,...
- Never place the electrode directly in the flask with the buffer fluid; it is best to pour a small amount into a separate, little container.
- Change electrodes:
 - if calibration is no longer possible;
 - if refilling (if applicable) is no longer possible;
 - if cleaning does not help;
 - if the electrode is damaged.

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